



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Rotor F-35-48-17, for 40 × conical 15 mL tubes

**€1,086.80**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Rotor F-35-48-17, for 40 × conical 15 mL tubes, incl. 48 steel sleeves and adapters for: Centrifuge 5810/5810 R, Centrifuge 5804/5804 R - High-capacity fixed-angle rotor for 15 mL conical tubes and glass tubes - Maximum speed: 5,005 × g (5,500 rpm)1) - Mold inserts for 15 mL conical tubes 1) Please note that older centrifuges require a software upgrade when using this rotor. For further information, please contact your local Eppendorf organization or your local Eppendorf distributor.

## Additional Information

---

No.	VB-2240
Manufacturer (Brand)	Eppendorf
VGKL number	106977200
EAN	4058072422752
Gross weight	4.6
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	4.40000 kg
Volume MAX	15000 µL
Country of origin	Germany
Packing width	0.32000 m
Packing height	0.22000 m
Packing depth	0.32000 m
Packaging volume	23000.000000 ccm
Customs tariff number	84219100
for model centrifuge	Eppendorf. 5804/5804 R
centrifugable up to	5005 g

