



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Rotor F-45-18-17-Cryo, without rotor lid, without adapter

**€516.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Rotor F-45-18-17-Cryo, without rotor lid, without adapter for: Centrifuge 5702/5702 R/5702 RH - Maximum speed: 1,970 × g (4,400 rpm) - Quiet centrifugation, even without rotor lid - Centrifugation of 1 - 2 mL cryo tubes (max. 13 mm Ø, use with adapter) and tubes with screw lid (max. 17 mm Ø, no adapter required) - Rotors and adapters are autoclavable (121 °C, 20 min)

## Additional Information

---

No.	VB-2045
Manufacturer (Brand)	Eppendorf
VGKL number	106947009
EAN	4058072420826
Gross weight	1.25
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	1.12500 kg
Volume MAX	2000 µl
Country of origin	Germany
Packing width	0.24000 m
Packing height	0.11000 m
Packing depth	0.24000 m
Packaging volume	7000.000000 ccm
Customs tariff number	84219100
centrifugable up to	1840 g

