




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor FA-10x5, aerosol-tight

€926.00
plus VAT &
Shipping

Product Images



Description

Rotor FA-10x5, aerosol-tight, incl. rotor lid, for 10 × 5.0 mL tubes (snap-on and screw-on lids) for: Centrifuge 5425 (obsolete), Centrifuge 5425 R (obsolete), Centrifuge 5425/5425 R - Max. Speed: 21,300 × g (15,060 rpm) - Max. Capacity: 10 × Eppendorf Tubes® 5.0 mL - Aerosol-tight with Eppendorf QuickLock® lid

Additional Information

No.	VB-1998
Manufacturer (Brand)	Eppendorf
VGKL number	106955008
EAN	4058072420369
Gross weight	1.145
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	1.03100 kg
Volume MAX	5000 µl
Country of origin	Germany
Customs tariff number	84219100
centrifugable up to	21330 g

