




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Rotor FA-18x2 kit, aerosol-tight

**€577.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Rotor FA-18x2 kit, aerosol-tight, incl. rotor lid, for 18 × 1.5/2.0 mL tubes or spin columns for: Centrifuge 5425 (obsolete), Centrifuge 5425 R (obsolete), Centrifuge 5425/5425 R - Max. Speed: 18,565 × g (15,060 rpm) - Rotor with raised rim for safe centrifugation of up to 18 spin columns - Aerosol-tight with Eppendorf QuickLock® lid

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | VB-2001          |
| Manufacturer (Brand)  | Eppendorf        |
| VGKL number           | 106978007        |
| EAN                   | 4058072420390    |
| Gross weight          | 1.258            |
| HSNumber              | 84219100         |
| Storage temperature   | Room temperature |
| Transport temperature | Room temperature |
| Net weight            | 1.13300 kg       |
| Volume MAX            | 2000 µl          |
| Country of origin     | Germany          |
| Customs tariff number | 84219100         |
| centrifugable up to   | 18565 g          |

