



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Rotor FA-20x5, incl. Eppendorf QuickLock® rotor lid

**€1,502.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

Rotor FA-20x5, aerosol-tight, for 20 Eppendorf Tubes® 5.0 mL, incl. Eppendorf QuickLock® rotor lid for: Centrifuge 5920 R, Centrifuge 5910 R, Centrifuge 5910 Ri – 20 places for Eppendorf Tubes® 5.0 mL1) – Max. Speed: 20,913 × g (13,100 rpm) – Eppendorf QuickLock rotor lid for fast lid locking – Aerosol-tight2) for safe centrifugation of hazardous samples – Anodized aluminum for high chemical resistance – Rotor, lid and adapter are autoclavable 1) For Eppendorf Tubes® 5.0 mL with screw cap, capacity is limited to 10 tubes. 2) Aerosol tightness tested by Public Health England, Porton Down (UK).

## Additional Information

---

No.	VB-2285
Manufacturer (Brand)	Eppendorf
EAN	4058072423186
Gross weight	3.2
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	2.88000 kg
Volume MAX	5000 µl
Country of origin	Germany
Customs tariff number	84219100
Fixed angle rotor	45°
Type rotor	Fixed angle rotor
for model centrifuge	Eppendorf. 5910 R
centrifugable up to	20913 g

