



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Rotor FA-24x2, aerosol-tight, incl. rotor lid, for 24 × 1.5/2.0 mL tubes

**€727.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Rotor FA-24x2, aerosol-tight, incl. rotor lid, for 24 × 1.5/2.0 mL tubes for: Centrifuge 5425 (obsolete), Centrifuge 5425/5425 R, Centrifuge 5425 R (obsolete) – Max. Speed: 21,300 × g (15,060 rpm) – Max. Capacity: 24 × 1.5/2.0 mL tubes – Aerosol-tight with Eppendorf QuickLock® lid

## Additional Information

---

No.	VB-1993
Manufacturer (Brand)	Eppendorf
VGKL number	106950006
EAN	4058072420338
Gross weight	1.23
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	1.10700 kg
Volume MAX	2000 µl
Country of origin	Germany
Customs tariff number	84219100
centrifugable up to	21330 g

