



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor FA-24x2, aerosol-tight, for 5420

€670.00
plus VAT &
Shipping

Product Images



Description

Rotor FA-24x2, aerosol-tight, incl. rotor lid, aerosol-tight, for 24 × 1.5/2.0 mL tubes, for 5420 for: Centrifuge 5420 - Max. Speed: 21,300 × g (15,060 rpm) - Max. Capacity: 24 × 1.5/2.0-mL tubes - Aerosol-tight

Additional Information

No.	VB-2533
Manufacturer (Brand)	Eppendorf
VGKL number	106940001
EAN	4058072425425
Gross weight	1.029
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.77200 kg
Country of origin	Germany
Customs tariff number	84219100
centrifugable up to	21300 g

