




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor FA-24x2-PTFE, aerosol-tight

€926.00
plus VAT &
Shipping

Product Images



Description

Rotor FA-24x2-PTFE, aerosol-tight, incl. rotor lid, PTFE-coated, with PTFE coating, for 24 × 1.5/2.0 mL tubes for: Centrifuge 5425 R (obsolete), Centrifuge 5425 (obsolete), Centrifuge 5425/5425 R – Max. Speed: 21,300 × g (15,060 rpm) – Max. Capacity: 24 × 1.5/2.0 mL tubes – Aerosol-tight with Eppendorf QuickLock® lid – PTFE coating for higher chemical resistance

Additional Information

No.	VB-1996
Manufacturer (Brand)	Eppendorf
VGKL number	106953005
EAN	4058072420345
Gross weight	1.235
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	1.11200 kg
Material	Polytetrafluoroethylene (PTFE)
Volume MAX	2000 µl
Country of origin	Germany
Customs tariff number	84219100
centrifugable up to	21330 g

