



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor FA-30x2, incl. aerosol-tight Eppendorf QuickLock® rotor lid

€1,068.00
plus VAT &
Shipping

Product Images



Description

Rotor FA-30x2, incl. aerosol-tight Eppendorf QuickLock® rotor lid for: Centrifuge 5910 Ri, Centrifuge 5910 R - 30 places for 1.5/2.0 mL tubes - Max. Speed: 20,984 × g (max. 13,700 rpm) - Eppendorf QuickLock® rotor lid for fast lid locking - Aerosol-tight for safe centrifugation of hazardous samples¹) - Anodized aluminum for high chemical resistance - Rotor and lid are autoclavable (121 °C, 20 min) ¹) Aerosol-tightness tested by Public Health England, Porton Down (UK).

Additional Information

No.	VB-2289
Manufacturer (Brand)	Eppendorf
EAN	4058072423223
Gross weight	2.2
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	2.05000 kg
Volume MAX	0.2 ml
Volume MIN	0.2 ml
Country of origin	Germany
Customs tariff number	84219100
for model centrifuge	Eppendorf. 5804/5804 R
centrifugable up to	20984 g

