



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor FA-45-20-17, incl. rotor lid, for 20 Eppendorf Tubes® 5.0 mL

€1,502.00
plus VAT &
Shipping

Product Images



Description

Rotor FA-45-20-17, incl. rotor lid, for 20 Eppendorf Tubes® 5.0 mL for: Centrifuge 5810/5810 R, Centrifuge 5804/5804 R – 20 places for Eppendorf Tubes® 5.0 mL) – Max. Speed: 20,913 × g (13,100 rpm)²⁾ – Aerosol-tight Eppendorf QuickLock® rotor lid³⁾ – Rotor, lid and adapter are autoclavable (121 °C, 20 min) – Anodized aluminum for high chemical resistance ¹⁾ For Eppendorf Tubes® 5.0 mL with screw cap, capacity is limited to 10 tubes. ²⁾ Please note that a software upgrade is required for older centrifuges when using this rotor. For more information, please contact your local Eppendorf organization or your local Eppendorf distributor. ³⁾ Aerosol tightness tested by Public Health England, Porton Down (UK).

Additional Information

No.	VB-2233
Manufacturer (Brand)	Eppendorf
EAN	4058072422684
Gross weight	3.2
HSNumber	84219100
Transport temperature	Room temperature
Storage temperature	Room temperature
Net weight	2.88000 kg
Volume MAX	5000 µl
Country of origin	Germany
Packing width	0.30000 m
Packing height	0.20000 m
Packing depth	0.30000 m
Packaging volume	18000.000000 ccm
Customs tariff number	84219100
Fixed angle rotor	45°
Type rotor	Fixed angle rotor
for model centrifuge	Eppendorf. 5804/5804 R
centrifugable up to	20913 g

