



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor FA-48x2, incl. aerosol-tight rotor lid Eppendorf QuickLock®

€1,502.00
plus VAT &
Shipping

Product Images



Description

Rotor FA-48x2, 48 × 1.5/2.0 mL tubes, incl. aerosol-tight rotor lid Eppendorf QuickLock® for: Centrifuge 5920 R, Centrifuge 5910 Ri, Centrifuge 5910 R – 48 places for 1.5/2.0 mL tubes – Max. Speed: 21,194 × g (max. 13,700 rpm) – Eppendorf QuickLock® rotor lid for fast lid closure – Aerosol-tight1) for safe centrifugation of hazardous samples – Anodized aluminum for high chemical resistance – Rotor and lid are autoclavable (121 °C, 20 min) 1) Aerosol-tightness tested by Public Health England, Porton Down (UK).

Additional Information

No.	VB-2286
Manufacturer (Brand)	Eppendorf
EAN	4058072423193
Gross weight	2.75
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	2.47500 kg
Volume MAX	2 ml
Volume MIN	0.2 ml
Country of origin	Germany
Customs tariff number	84219100
for model centrifuge	Eppendorf. 5804/5804 R
centrifugable up to	21194 g

