



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Rotor FA-6x50, aerosol-tight, incl. rotor lid, 6 × 50 mL conical tubes

€2,448.00
plus VAT &
Shipping

Product Images



Description

Rotor FA-6x50, aerosol-tight, incl. rotor lid, 6 × 50 mL conical tubes for: Centrifuge 5910 R, Centrifuge 5910 Ri, Centrifuge 5920 R - High speed rotor for applications with tubes from 5 - 50 mL - Max. Speed 20,130 × g (12,100 rpm) - Eppendorf QuickLock® rotor lid for fast lid locking - Aerosol-tight¹) for safe centrifugation of hazardous samples - Low rotor weight (3.450 kg) for easy handling - Anodized aluminum for high chemical resistance - 45° bore angle minimizes pellet streaks - Rotor, lid and adapter are autoclavable ¹) Aerosol-tightness tested by Public Health England, Porton Down (UK).

Additional Information

No.	VB-2287
Manufacturer (Brand)	Eppendorf
EAN	4058072423209
Gross weight	4.05
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	3.50000 kg
Country of origin	Germany
Customs tariff number	84219100
Fixed angle rotor	45°
for model centrifuge	Eppendorf. 5910 R

