



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Rotor T-60-11, incl. rotor lid, for 1.5/2.0 mL tubes, without adapter

**€1,740.20**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Rotor T-60-11, incl. rotor lid, for 1.5/2.0 mL tubes, without adapter for: Centrifuge 5810/5810 R, Centrifuge 5804/5804 R - High capacity for 60 × 1.5/2.0-mL tubes - Max. Speed: 16,435 × g (14,000 rpm) - Horizontal centrifugation for pellet formation at the bottom of the tube - Ideal for small sample volumes

## Additional Information

---

No.	VB-2079
Manufacturer (Brand)	Eppendorf
VGKL number	556930003
EAN	4058072421168
Gross weight	2.9
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	2.61000 kg
Country of origin	Germany
Packing width	0.24000 m
Packing height	0.15000 m
Packing depth	0.24000 m
Packaging volume	9000.000000 ccm
Customs tariff number	84219100
Fixed angle rotor	90°
for model centrifuge	Eppendorf. 5804/5804 R

