



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Rotor f. Pico + Fresco 21: 24 x 1.5 ml, with aerosol-tight lid

€719.00
plus VAT &
Shipping

Product Images



Description

Order No. 3-1885

Aerosol-tight ClickSeal™ rotor offers first-class safety and convenience
Rotor for 24 x 1.5/2.0 ml with aerosol-tight ClickSeal biocontainment lid

Additional Information

No.	LD-1998
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1885
EAN	4058072287733
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Rotor(s) is/are suitable for	Reaction tubes. 24 x 1.5/ 2.0 ml
Type rotor	Fixed angle rotor
for model centrifuge	Thermo Scientific Heraeus Fresco 17 Thermo Scientific Heraeus Fresco 21 Thermo Scientific Heraeus Pico 21 Thermo Scientific Heraeus Pico 17

