



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Fixed Angle Rotor f. 10x15 ml Falcon, incl. 10 adapters & 10 rubber mats for E-1902

**€517.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Fixed-angle rotors, Eppendorf

Max. Speed: 4,400 rpm, autoclavable.

#### **Order no. E-1906**

Fixed-angle rotor type F-35-30-17

for 10 x 15 ml tubes, incl. 10 adapters and 10 rubber mats

## Additional Information

---

No.	VB-2016
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-1906
EAN	4058072420543
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Rotor(s) is/are suitable for	Centrifuge tubes. 10 x15 ml
Type rotor	Fixed angle rotor
for model centrifuge	Eppendorf. 5702 Eppendorf. 5702 R Eppendorf. 5702 RH

