



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Fixed Angle Rotor f. 30 x 1.5/2.0 ml reaction tubes, incl. aluminum lid

**€729.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Fixed-angle rotors, Eppendorf

autoclavable

#### **Order No. E-1896**

Fixed-angle rotor type F-45-30-11

for 30 x 1.5/2.0 ml reaction tubes, incl. aluminum lid, max. 14,000 rpm

## Additional Information

---

No.	VB-2069
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-1896
EAN	4058072421069
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Rotor(s) is/are suitable for	Reaction tubes. 30 x 1.5/2.0 ml
Type rotor	Fixed angle rotor
for model centrifuge	Eppendorf. 5804 Eppendorf. 5804 R

