



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® safety net red, for DM 50-100 mm, 5 meters

**€40.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

neoLab **safety nets made of PE in** various colors and sizes.

Ideally suited as shatter protection for covering various glass containers in the laboratory, such as bottles, flasks, desiccators and other vacuum glassware. This reduces the range of flying glass fragments in case of explosion or implosion.

Easy to cut using scissors. 5 m / roll.

## Additional Information

---

No.	2-3094
Manufacturer (Brand)	neoLab
EAN	4058072112417
Transport temperature	Room temperature
Color	red
Material	Polyethylene (PE)
Wide	100 mm
Length	5 m
suitable for	Bottles. round Bottles. square
TBST MAX	80 °C
TBST MIN	-50 °C

