



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® safety net yellow, for DM 90-200 mm, 5 meters

€54.00
plus VAT &
Shipping

Product Images



Description

neoLab **safety nets made of PE in** various colors and sizes.

Ideally suited as shatter protection for covering various glass containers in the laboratory, such as bottles, flasks, desiccators and other vacuum glassware. This reduces the range of flying glass fragments in case of explosion or implosion.

Easy to cut using scissors. 5 m / roll.

Additional Information

No.	2-3095
Manufacturer (Brand)	neoLab
EAN	4058072112424
Transport temperature	Room temperature
Color	yellow
Material	Polyethylene (PE)
Wide	200 mm
Length	5 m
suitable for	Bottles. round Bottles. square
TBST MAX	80 °C
TBST MIN	-50 °C

