



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN RS-232 interface cable for connecting an external device

€50.00
plus VAT &
Shipping

Product Images



Description

- Cable length approx. 1.5 m

Suitable for

- Precision balances KERN 440, KERN 572, KERN 573, KERN KB, KERN PCB, KERN PCD, PLS-A, PLJ-A, KERN PLT
- Analytical balances ALS-A, KERN ALT-B
- Floor scales KERN BTBP, KERN BTSP, KERN BTEP
- Moisture analyzers KERN DBS, KERN DLB
- Counting scales KERN CDS, KERN CDE, KERN CKE
- Platform scales KERN DE, KERN DS, KERN IKT, KERN ITB, KERN ITS, KERN ITT, KERN NDE
- Bench scales KERN FCB, KERN FTB, KERN FTC, KERN FKB-A, KERN FKB, KERN FKT

Additional Information

No.	KP-0518
Manufacturer (Brand)	KERN
EAN	4058072125936
Transport temperature	Room temperature
Manufacturer	CORE

