



neolabLine pH test strips

Product Images



Description

The neoLab Universal **pH indicator strips** are ideally suited for a quick and meaningful pH value determination from pH 0-14. Four test fields on the strip in conjunction with a color comparison scale enable an exact determination of the pH value. A basis for decision-making is created within a few seconds. Packed 100 to a handy plastic square pack, the neoLab pH indicator sticks can be taken anywhere and used outside the laboratory.

- Specifications:
- pH 0-14, Graduation: 1.0
 - No bleeding
 - Good readability
 - No contamination of the medium

Additional Information

No.	3-0000
Manufacturer (Brand)	neoLabLine
EAN	4058072115579
Wide	6 mm
Length	80 mm
pH MAX	14.0 pH
pH MIN	0.0 pH
Format/ Packing	einzeln. im Karton/ Spender
Sensitivity	1 pH
pH levels	14

