



## Macherey & Nagel pH-Fix 6.0-7.7

### Product Images



### Description

pH-Fix 6.0-7.7 Indicator strips Measuring range: pH 6.0-6.4-6.7-7.0-7.3-7.7 Pg. of 100 test strips 5.5 x 85 mm

## Additional Information

|                       |                     |
|-----------------------|---------------------|
| No.                   | ML-6717             |
| Dual use indicator    | false               |
| Manufacturer (Brand)  | Macherey & Nagel    |
| VGKL number           | 244692150           |
| EAN                   | 4058072175498       |
| Gross weight          | 0.0319              |
| Manufacturer EAN      | 4046681001512       |
| Storage temperature   | Raumtemperatur      |
| Net weight            | 0.03190 kg          |
| S-phrases             | S 22                |
| Country of origin     | Deutschland         |
| Region of origin      | Nordrhein-Westfalen |
| Packing width         | 0.06500 m           |
| sterile               | No                  |
| Packing height        | 0.09000 m           |
| Packing depth         | 0.01000 m           |
| Packaging volume      | 0.00005850 ccm      |
| Customs tariff number | 38229000            |
| pH MAX                | 7.7 pH              |
| pH MIN                | 6.0 pH              |
| Format/ Packing       | Stäbchen            |

