



Macherey & Nagel pH-Fix 7.0-14.0

Product Images



Description

pH-Fix 7.0-14.0 Indicator strips Measuring range: pH 7.0-7.5-8.0-8.5-9.0-9.5-10.0- 10.5-11.0-11.5-12.0-12.5-13.0-13.5-14.0 Pg. à 100 test strips 5.5 x 85 mm

Additional Information

No.	ML-6711
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	244692125
EAN	4058072166014
Gross weight	0.0345
Manufacturer EAN	4046681001451
Storage temperature	Raumtemperatur
Net weight	0.03450 kg
S-phrases	S 22
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.06600 m
sterile	No
Packing height	0.09000 m
Packing depth	0.01000 m
Packaging volume	0.00005900 ccm
Customs tariff number	38229000
pH MAX	14.0 pH
pH MIN	7.0 pH
Format/ Packing	Stäbchen

