



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Hot Air Sterilization Indicators, 1000 pcs./pack

€173.00
plus VAT &
Shipping

Product Images



Description

Indicators for sterilization with hot air, e.g. syringes, flasks, bottles, pipette tubes or dishes.

Color change from green to black at +160°C to +180°C.

Items to be stored should be wrapped in aluminum foil beforehand and a second label applied on the inside.

Note: Store at normal room temperature of 10° - 38°C (50° - 100°F) and 10-70% relative humidity. No special storage conditions are required after exposure to dry heat

Additional Information

No.	3-7410
Manufacturer (Brand)	neoLab
EAN	4058072020170
Transport temperature	Room temperature
suitable for	Hot air sterilization
Format/ Packing	Labels

