



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Steri control tubes for autoclaves up to 125°C , 100 pcs./pack

€117.60
plus VAT &
Shipping

Product Images



Description

Correct sterilization is essential, especially in medical and bacteriological laboratories, but is not always guaranteed. The color change of the glass control tube from red to green indicates correct sterilization. Place the tubes as close as possible to or in the load. Store in a cool and dry place! Suitable for autoclaves up to 125°C.

Shelf life: 6 months

Additional Information

No.	3-1280
Manufacturer (Brand)	neoLab
EAN	4058072012731
Transport temperature	Room temperature
diameter	6 mm
Length	45 mm
suitable for	Steam sterilization (in autoclave)

