




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Rasotherm beaker high form with spout, (Boro 3.3), 50 ml

Special Price

€1.90 was

€2.47 As low as

€1.86

**plus VAT &
Shipping**

Product Images



Description

The RASOTHERM® **beakers**, high form with spout are made of especially heat & chemical resistant borosilicate glass 3.3.

With graduation.

All beakers meet the highest quality standards and are manufactured according to the specifications of the current ISO standard 3819.

Thus, the beakers are suitable for many purposes in your laboratory.

Additional Information

| | |
|----------------------------|--|
| No. | 1-0164 |
| Manufacturer (Brand) | Rasotherm |
| EAN | 4058072080983 |
| autoclavable at 121°C | Yes |
| diameter | 38 mm |
| Color | transparent (clear) |
| sterile | No |
| Form | high shape |
| Height | 70 mm |
| Cups and scales properties | with spout graduated suitable for foodstuffs |
| Material | Borosilicate glass |
| Color graduation | white |
| Material cup | Borosilicate glass |
| Type cup | Beakers |
| Volume cup | 50 ml |

