



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Rasotherm beaker high form with spout, (Boro 3.3), 250 ml

Special Price

**€2.30** was


**€2.97** As low as

**€2.23**

**plus VAT &  
Shipping**

## Product Images

---



## Description

---

The RASOTHERM® **beakers**, high form with spout are made of especially heat & chemical resistant borosilicate glass 3.3.

With graduation.

All beakers meet the highest quality standards and are manufactured according to the specifications of the current ISO standard 3819.

Thus, the beakers are suitable for many purposes in your laboratory.

## Additional Information

---

No.	1-0167
Manufacturer (Brand)	Rasotherm
EAN	4058072080891
autoclavable at 121°C	Yes
diameter	60 mm
Color	transparent (clear)
sterile	No
Form	high shape
Height	120 mm
Cups and scales properties	with spout graduated suitable for foodstuffs
Material	Borosilicate glass
Color graduation	white
Material cup	Borosilicate glass
Type cup	Beakers
Volume cup	250 ml

