



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Erlenmeyer flask, graduated, GL 25, with cap (PBT), 100 ml

€400.60
plus VAT &
Shipping

Product Images



Description

DURAN® Erlenmeyer flasks, with DIN thread, with PBT closure, 100 ml The original DURAN® Erlenmeyer flasks are predestined for use in laboratories due to their excellent properties, such as very good chemical resistance, virtually inert behavior and a high service temperature. In addition, due to their uniform wall thickness distribution, they are characterized by minimal thermal expansion and the resulting high thermal shock resistance. DURAN® Erlenmeyer flasks are therefore particularly suitable for work at very high temperatures and rapid temperature changes. Special properties and advantages: The conical shape facilitates mixing Easy-to-read scale with graduation for measuring the contents Baked-on and therefore permanently durable white printed image The flask can be closed with a PBT closure or membrane closure (gas exchange possible) Application: Cultivation, storage and transport, mixing of several components, heating, diluting solutions, measuring liquid quantities, supplied with or without closure.

Additional Information

No.	SD-0826
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	102115324
EAN	4058072345570
Manufacturer EAN	4032051014954
Gross weight	0.084
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.07000 kg
Material	Borosilicate glass 3.3
Volume	100 ml
Country of origin	Germany
sterile	No
Region of origin	Rhineland-Palatinate
Packaging volume	0.00469500 cc
Customs tariff number	70172000
Thread	GL 25

