



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Rasotherm Erlenmeyer flask wide neck, (Boro 3.3), 200 ml

Special Price

**€4.50** was

**€5.82** As low as

**€4.08**

**plus VAT &  
Shipping**

## Product Images

---



## Description

---

The RASOTHERM® **Erlenmeyer flasks** with wide neck are made of particularly heat- & chemical-resistant borosilicate glass 3.3.

All flasks meet the highest quality standards and comply with ISO 24450.

This makes Erlenmeyer flasks suitable for many applications in your laboratory.

All flasks are graduated.

## Additional Information

---

No.	1-0188
Manufacturer (Brand)	Rasotherm
EAN	4058072081225
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	200 ml
sterile	No
Height	130 mm
Type closure/ lid	without closure
Pistons properties	with certification graduated
Color graduation	white
Neck size	far
DM Neck	60 mm
DM Piston	79 mm
Type piston	Erlenmeyer flask
Type edge	Flared rim
Type certificate	Declaration of conformity (standardization)
standardized according to ISO standard	ISO 1773

