



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Rasotherm Erlenmeyer flask wide neck, (Boro 3.3), 300 ml

Special Price

€4.20 was

€5.42 As low as

€3.83

**plus VAT &
Shipping**

Product Images



Description

The RASOTHERM® **Erlenmeyer flasks** with wide neck are made of particularly heat- & chemical-resistant borosilicate glass 3.3.

All flasks meet the highest quality standards and comply with ISO 24450.

This makes Erlenmeyer flasks suitable for many applications in your laboratory.

All flasks are graduated.

Additional Information

No.	1-0190
Manufacturer (Brand)	Rasotherm
EAN	4058072081126
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	300 ml
sterile	No
Height	156 mm
Type closure/ lid	without closure
Pistons properties	with certification graduated
Color graduation	white
Neck size	far
DM Neck	60 mm
DM Piston	91 mm
Type piston	Erlenmeyer flask
Type edge	Flared rim
Type certificate	Declaration of conformity (standardization)
standardized according to ISO standard	ISO 1773

