



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Rasotherm Erlenmeyer flask wide neck, (Boro 3.3), 2000 ml

Special Price

€13.90 was

€17.80 As low as

€12.90

**plus VAT &
Shipping**



Description

The RASOTHERM® **Erlenmeyer flasks** with wide neck are made of particularly heat- & chemical-resistant borosilicate glass 3.3.

All flasks meet the highest quality standards and comply with ISO 24450.

This makes Erlenmeyer flasks suitable for many applications in your laboratory.

All flasks are graduated.

Additional Information

| | |
|--|---|
| No. | 1-0193 |
| Manufacturer (Brand) | Rasotherm |
| EAN | 4058072081270 |
| autoclavable at 121°C | Yes |
| Color | transparent (clear) |
| Material | Borosilicate glass |
| Volume | 2000 ml |
| sterile | No |
| Height | 310 mm |
| Type closure/ lid | without closure |
| Pistons properties | with certification graduated |
| Color graduation | white |
| Neck size | far |
| DM Neck | 69 mm |
| DM Piston | 187 mm |
| Type piston | Erlenmeyer flask |
| Type edge | Flared rim |
| Type certificate | Declaration of conformity (standardization) |
| standardized according to ISO standard | ISO 1773 |

