



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Rasotherm Erlenmeyer flask narrow neck, (Boro 3.3), 1000 ml

Special Price

€6.60 was

€8.40 As low as

€5.94

**plus VAT &
Shipping**

Product Images



Description

The RASOTHERM® **Erlenmeyer flasks** with narrow neck are made of particularly heat- & chemical-resistant borosilicate glass 3.3.

With graduation.

All flasks meet the highest quality standards and are manufactured according to the current ISO standard 1773. This makes the Erlenmeyer flasks suitable for many applications in your laboratory.

Additional Information

No.	1-0181
Manufacturer (Brand)	Rasotherm
EAN	4058072081157
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	1000 ml
sterile	No
Height	220 mm
Pistons properties	graduated
DM Neck	38 mm
DM Piston	131 mm
Color graduation	white
Neck size	narrow
Type piston	Erlenmeyer flask
Type edge	Flared rim
standardized according to ISO standard	ISO 1773

