



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Rasotherm Erlenmeyer flask narrow neck, (Boro 3.3), 3000 ml


Special Price

**€27.50** was

**€35.20**

**plus VAT &  
Shipping**

### Product Images



## Description

---

The RASOTHERM® **Erlenmeyer flasks** with narrow neck are made of particularly heat- & chemical-resistant borosilicate glass 3.3.

With graduation.

All flasks meet the highest quality standards and are manufactured according to the current ISO standard 1773. This makes the Erlenmeyer flasks suitable for many applications in your laboratory.

## Additional Information

---

No.	1-0183
Manufacturer (Brand)	Rasotherm
EAN	4058072081171
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	3000 ml
sterile	No
Height	310 mm
Pistons properties	graduated
DM Neck	50 mm
DM Piston	187 mm
Color graduation	white
Neck size	narrow
Type piston	Erlenmeyer flask
Type edge	Flared rim
standardized according to ISO standard	ISO 1773

