



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz iodine number flask with standard ground joint, without collar, 250 ml, with hexagonal hollow stopper

**€190.10
plus VAT &
Shipping**

Product Images



Description

Iodine number flask with standard ground joint, without collar.

With graduation for reading the approximate volume and hexagonal hollow stopper NS 29/32.

With labeling field.

Additional Information

No.	LZ-0525
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2525
EAN	4058072302672
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass DURAN
Volume	250 ml
sterile	No
Height	160 mm
Type closure/ lid	NS lid plug
Pistons properties	with closure/lid with standard ground graduated
DM Piston	85 mm
Color graduation	white
Standard ground	29/32
Type piston	Iodine number flask
Type edge	Standard ground

