



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand volumetric flask cl. A blue, w. NS-PP
stopper, 500 ml, 2 pcs./pack**

**€61.60
plus VAT &
Shipping**

Product Images



Description

BLAUBRAND®, conformity certified.

DURAN® DIN EN ISO 1042. Adjusted to "In". Ring mark and labeling in high-contrast blue enamel color.

All BLAUBRAND® volumetric flasks are supplied with batch number and an enclosed batch certificate per original packaging.

With the exception of the wide-neck version, the volumetric flasks are also available with USP individual and batch certificate.

Additional Information

No.	BN-1115
Manufacturer (Brand)	Brand
old neoLab article no.	4-4068
EAN	4058072205164
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Material	Borosilicate glass
Volume	500 ml
Form	around
Type closure/ lid	NS plug
Pistons properties	with certification with closure/lid with ring mark with standard ground joint
Color lettering	blue
Color ring mark	blue
Class	A
Material closure/ lid	Polypropylene (PP)
Standard ground	19/26
Type piston	Volumetric flask
Type certificate	Batch certificate from the manufacturer
standardized according to ISO standard	ISO 1042

