



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand volumetric flask cl. A blue, w. NS-PP  
stopper, 1000 ml, 2 pcs./pack**

**€54.40  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

BLAUBRAND®, conformity certified.

DURAN® DIN EN ISO 1042. Adjusted to "In". Ring mark and labeling in high-contrast blue enamel color.

All BLAUBRAND® volumetric flasks are supplied with batch number and an enclosed batch certificate per original packaging.

With the exception of the wide-neck version, the volumetric flasks are also available with USP individual and batch certificate.

## Additional Information

No.	BN-1116
Manufacturer (Brand)	Brand
old neoLab article no.	4-4069
EAN	4058072205171
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Material	Borosilicate glass
Volume	1000 ml
Form	around
Type closure/ lid	NS plug
Pistons properties	with certification with closure/lid with ring mark with standard ground joint
Class	A
Standard ground	24/29
Color lettering	blue
Color ring mark	blue
Material closure/ lid	Polypropylene (PP)
Type piston	Volumetric flask
Type certificate	Batch certificate from the manufacturer
standardized according to ISO standard	ISO 1042

