



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hirschmann® Volumetric flasks Duran glass, PE stopper, NS 14/23, 100 ml, 2 pcs./pack

€46.70
plus VAT &
Shipping

Product Images



Description

Volumetric flask with NS plastic stopper
Class A, Duran glass
2 pieces per pack

DIN 12 664

Additional Information

No.	NH-0371
Manufacturer (Brand)	Hirschmann
old neoLab article no.	E-1569
EAN	4058072336363
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Material	Borosilicate glass
Volume	100 ml
Form	round
Type closure/ lid	NS plug
Pistons properties	with certification with closure/lid with ring mark with standard ground joint
Class	A
Standard ground	14/23
Color lettering	blue
Color ring mark	blue
Material closure/ lid	Low-density polyethylene (LDPE)
Type piston	Volumetric flask
Type certificate	Batch certificate from the manufacturer
standardized according to ISO standard	ISO 1042

