



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® volumetric flask trapezoidal 10 ml, NS stopper made of PE 7/16, 2 pcs./pack

€33.60
plus VAT &
Shipping

Product Images



Description

Adjusted to IN. Class A. With certificate of conformity. With standard ground joint (NS 7/16) and poly stopper. Ring mark in brown color. With wide, flat base for optimum stability, DIN EN ISO 1042. DE-M marking.

Made of borosilicate glass 3.3

Additional Information

No.	E-1692
Manufacturer (Brand)	neolab
EAN	4058072066574
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Material	Borosilicate glass
Volume	10 ml
Form	trapezoidal
Type closure/ lid	NS plug
Pistons properties	with certification with closure/lid with ring mark with standard ground joint
Class	A
Standard ground	7/16
Color lettering	brown
Color ring mark	brown
Material closure/ lid	High-density polyethylene (HDPE)
Type piston	Volumetric flask
Type certificate	Batch certificate from the manufacturer
standardized according to ISO standard	ISO 1042

