



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® volumetric flask trapezoidal 20 ml, NS stopper made of PE 10/19, 2 pcs./pack

**€33.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Adjusted to IN. Class A. With certificate of conformity. With standard ground joint (NS 10/19) and poly stopper. Ring mark in brown color. With wide, flat base for optimum stability, DIN EN ISO 1042. DE-M marking.

Made of borosilicate glass 3.3

## Additional Information

---

No.	E-1693
Manufacturer (Brand)	neolab
EAN	4058072066581
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Material	Borosilicate glass
Volume	20 ml
Form	trapezoidal
Type closure/ lid	NS plug
Pistons properties	with certification with closure/lid with ring mark with standard ground joint
Color lettering	brown
Color ring mark	brown
Class	A
Material closure/ lid	High-density polyethylene (HDPE)
Standard ground	10/19
Type piston	Volumetric flask
Type certificate	Batch certificate from the manufacturer
standardized according to ISO standard	ISO 1042

