



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® volumetric flask trapezoidal 25 ml, NS
stopper made of PE 10/19, 2 pcs./pack**

**€33.70
plus VAT &
Shipping**

Product Images



Description

Adjusted to IN. Class A. With certificate of conformity. With standard ground joint (NS 10/19) and poly stopper. Ring mark in brown color. With wide, flat bottom for optimum stability, DIN EN ISO 1042. DE-M marking.
Made of borosilicate glass 3.3

Additional Information

No.	E-1694
Manufacturer (Brand)	neoLab
EAN	4058072066598
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Material	Borosilicate glass
Volume	25 ml
Form	trapezoidal
Type closure/ lid	NS plug
Pistons properties	with certification with closure/lid with ring mark with standard ground joint
Class	A
Standard ground	10/19
Color lettering	brown
Color ring mark	brown
Material closure/ lid	High-density polyethylene (HDPE)
Type piston	Volumetric flask
Type certificate	Batch certificate from the manufacturer
standardized according to ISO standard	ISO 1042

