



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




## Brand measuring cylinder high form, TPX, cl. A, blue grade, 50 ml

**€12.87** As low as  
**€9.90**  
plus VAT &  
Shipping

## Product Images

---



## Description

---

PLASTIBRAND® PMP is crystal clear and has excellent chemical resistance.

Blue printed graduation, extremely durable.

Increased stability due to hexagonal base with knobs.

Extremely low taper allows even graduation.

Functional spout.

Thermal loads up to +121°C (autoclaving) do not permanently exceed the error limit.

To protect the graduation, cleaning up to max. +60°C is recommended.

### **Order No. 4-4112**

Class A: according to DIN12681 and ISO 6706, with batch certificate

## Additional Information

No.	BN-0993
Manufacturer (Brand)	Brand
old neoLab article no.	4-4112
EAN	4058072204013
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polymethylpentene (PMP. TPX)
Volume	50 ml
Form	high shape
Height	200 mm
Color graduation	blue
Class	A
Division	1 ml
DM Neck	29 mm
Type certificate	Declaration of conformity (standardization)
standardized according to DIN standard (DIN)	DIN 12681
standardized according to ISO standard	ISO 6706

