



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657




Brand measuring cylinder high form, TPX, cl. A, blue grade, 250 ml

€27.17 As low as

€20.90

**plus VAT &
Shipping**

Product Images



Description

PLASTIBRAND® PMP is crystal clear and has excellent chemical resistance.

Blue printed graduation, extremely durable.

Increased stability due to hexagonal base with knobs.

Extremely low taper allows even graduation.

Functional spout.

Thermal loads up to +121°C (autoclaving) do not permanently exceed the error limit.

To protect the graduation, cleaning up to max. +60°C is recommended.

Order No. 4-4114

Class A: according to DIN 12681 and ISO 6706, with batch certificate

Additional Information

No.	BN-0995
Manufacturer (Brand)	Brand
old neoLab article no.	4-4114
EAN	4058072204037
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polymethylpentene (PMP. TPX)
Volume	250 ml
Form	high shape
Height	315 mm
DM Neck	53 mm
Color graduation	blue
Class	A
Division	2 ml
Type certificate	Declaration of conformity (standardization)
standardized according to DIN standard (DIN)	DIN 12681
standardized according to ISO standard	ISO 6706

