





neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand measuring cylinder high form, TPX, cl. A, blue grade, 500 ml

€41.21 As low as
€31.70
plus VAT &
Shipping



Description

PLASTIBRAND® PMP is crystal clear and has excellent chemical resistance.

Blue printed graduation, extremely durable.

Increased stability due to hexagonal base with knobs.

Extremely low taper allows even graduation.

Functional spout.

Thermal loads up to +121°C (autoclaving) do not permanently exceed the error limit.

To protect the graduation, cleaning up to max. +60°C is recommended.

Order No. 4-4115

Class A: according to DIN 12681 and ISO 6706, with batch certificate

Additional Information

No.	BN-0996
Manufacturer (Brand)	Brand
old neoLab article no.	4-4115
EAN	4058072204044
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polymethylpentene (PMP. TPX)
Volume	500 ml
Form	high shape
Height	360 mm
Color graduation	blue
Class	A
Division	5 ml
DM Neck	65 mm
Type certificate	Declaration of conformity (standardization)
standardized according to DIN standard (DIN)	DIN 12681
standardized according to ISO standard	ISO 6706

