



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand measuring cylinder high form, TPX, cl. B, blue grade, 50 ml

€10.68 As low as
€8.90
plus VAT &
Shipping

Product Images



Description

PLASTIBRAND® PMP is crystal clear and has excellent chemical resistance.
Blue printed graduation, extremely durable.
Increased stability due to hexagonal base with knobs.
Extremely low taper allows even graduation.
Functional spout.
Thermal loads up to +121°C (autoclaving) do not permanently exceed the error limit.
To protect the graduation, cleaning up to max. +60°C is recommended.

Order No. 4-4102

Class B: according to DIN 12681 and ISO 6706

Additional Information

No.	BN-0961
Manufacturer (Brand)	Brand
old neoLab article no.	4-4102
EAN	4058072203771
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polymethylpentene (PMP. TPX)
Volume	50 ml
Form	high shape
Height	200 mm
Color graduation	blue
Class	B
Division	1 ml
DM Neck	29 mm
standardized according to DIN standard (DIN)	DIN 12681
standardized according to ISO standard	ISO 6706

