



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




**Brand measuring cylinder 10 ml DURAN glass,  
high form, cl. B, silver Eterna brown, 2 pcs./pack**

**€19.30  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

SILVER FIRING ETERNA, DURAN®. DIN EN ISO 4788.

Adjusted to "In". With spout and hexagonal base. Graduation. Ring mark and lettering in brown ETERNA diffusion color.

## Additional Information

---

No.	BN-0910
Manufacturer (Brand)	Brand
old neolab article no.	4-4051
EAN	4058072186500
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	10 ml
Form	high shape
Height	140 mm
DM Neck	14 mm
Color graduation	brown
Class	B
Division	0.2 ml
standardized according to DIN standard (DIN EN ISO)	DIN EN ISO 4788:2005-08

