



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand measuring cylinder 25 ml DURAN glass,
high form, cl. B, silver Eterna brown, 2 pcs./pack**

**€20.20
plus VAT &
Shipping**

Product Images



Description

SILVER FIRING ETERNA, DURAN®. DIN EN ISO 4788.

Adjusted to "In". With spout and hexagonal base. Graduation. Ring mark and lettering in brown ETERNA diffusion color.

Additional Information

No.	BN-0911
Manufacturer (Brand)	Brand
old neolab article no.	4-4052
EAN	4058072168070
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	25 ml
Form	high shape
Height	170 mm
Color graduation	brown
Class	B
Division	0.5 ml
DM Neck	21 mm
standardized according to DIN standard (DIN EN ISO)	DIN EN ISO 4788:2005-08

