



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




---

## Brand measuring cylinder 2000 ml DURAN glass, high form, cl. B, silver Eterna brown

**€100.05**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

SILVER FIRING ETERNA, DURAN®. DIN EN ISO 4788.

Adjusted to "In". With spout and hexagonal base. Graduation. Ring mark and lettering in brown ETERNA diffusion color.

## Additional Information

---

No.	BN-0917
Manufacturer (Brand)	Brand
old neolab article no.	4-4058
EAN	4058072203412
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	2000 ml
Form	high shape
Height	505 mm
Color graduation	brown
Class	B
Division	20 ml
DM Neck	85 mm
standardized according to DIN standard (DIN EN ISO)	DIN EN ISO 4788:2005-08

