



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Brand graduated cylinder cl. B, 250 ml, low  
form, silver Eterna brown, 2 pcs/pack**

**€44.40  
plus VAT &  
Shipping**



## Description

---

SILVER FIRING ETERNA, DURAN®. DIN EN ISO 4788. adjusted to "In".

With spout and hexagonal base. Graduation. Ring mark and lettering in brown ETERNA diffusion color.

## Additional Information

---

No.	BN-1227
Manufacturer (Brand)	Brand
old neolab article no.	4-4045
EAN	4058072206079
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	250 ml
Form	low form
Height	195 mm
Color graduation	brown
Class	B
Division	5 ml
DM Neck	54 mm
standardized according to DIN standard (DIN EN ISO)	DIN EN ISO 4788:2005-08

