




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand graduated cylinder cl. B, 500 ml, low form, silver Eterna brown, 2 pcs/pack

**€69.00
plus VAT &
Shipping**



Description

SILVER FIRING ETERNA, DURAN®. DIN EN ISO 4788. adjusted to "In".

With spout and hexagonal base. Graduation. Ring mark and lettering in brown ETERNA diffusion color.

Additional Information

No.	BN-1228
Manufacturer (Brand)	Brand
old neoLab article no.	4-4046
EAN	4058072206086
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
Form	low form
Height	250 mm
Color graduation	brown
Class	B
Division	10 ml
DM Neck	65 mm
standardized according to DIN standard (DIN EN ISO)	DIN EN ISO 4788:2005-08

