



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand mixing cylinder, high form, SILBERBRAND
ETERNA, class B, Boro 3.3, 25 ml: 0.5ml, NS 14/23,
PP stopper**

**€40.70
plus VAT &
Shipping**



Description

Boro 3.3. DIN EN ISO 4788. Adjusted to "In". With PP stopper and hexagonal base. Height indication without stopper.

Additional Information

No.	BN-0953
Manufacturer (Brand)	Brand
EAN	4058072175832
Gross weight	0.1704
HSNumber	70172000
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.13140 kg
Material	Borosilicate glass 3.3
Volume	25 ml
Country of origin	Germany
Packing width	0.06800 m
Packing height	0.05700 m
Packing depth	0.25200 m
Packaging volume	976.75200000 cc
Customs tariff number	70172000
Form	high shape
Color graduation	silver
Class	B
Standard ground	14/23
Division	0.5 ml
adjusted to	"in"

