



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz receiving flask with KS bowl S 35, 2,000 ml

€91.20
plus VAT &
Shipping

Product Images



Description

Especially for rotary evaporators, with spherical ground joint bowl S 35, with labeling field, material: DURAN®

Additional Information

No.	LZ-0578
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2578
EAN	4058072303181
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	166 mm
Color	transparent (clear)
Volume	2000 ml
Form	spherical
Pistons properties	with cut
Spherical cut bowl	S35
Type glass jar	Receiving flask for rotary evaporator
Type ground	Spherical ground

